



DISPOSABLE MATRIX
SYSTEM

Frequently Asked Questions

What is the difference between the four colors?

Green: Pro-Matrix Narrow (4.5mm) for pre-molars and children

Blue: Pro-Matrix Wide (6mm) for molars

Yellow: Pro-Matrix Curve Narrow (4.5mm) for pre-molars and children

Orange: Pro-Matrix Curve Wide (6mm) for molars

What is the advantage of Pro-Matrix Curve?

The bands are pre-contoured for more natural looking restorations and more predictable contact points.

Is it delivered sterile?

No. It is clean, equivalent to other types of matrix systems. Each Pro-Matrix is individually packed.

How thick is the Pro-Matrix band?

0.04mm

Can it be used with wedges?

Yes. Any wedge can be used. See image below.



Can it be used with both composite and amalgam?

Yes. Pro-Matrix is suitable for both due to its high-tension mechanism allowing for compacting of amalgam.

Will it burnish?

No. This product is a single-use disposable item.

Can it be autoclaved?

Yes. Pro-Matrix is suitable for both due to its high-tension mechanism allowing for compacting of amalgam.

Can more than one Pro-Matrix be used next to each other?

Yes. They stack tightly for patient comfort and excellent visibility.




THE ORIGINAL

**TRUSTED IN 20,000,000+
PROCEDURES & COUNTING!**

For more information, email customerservice@younginnovations.com or call (800) 325-1881.

Instructions for Use



The toggle should be positioned close to the gum line to ensure the correct fit around the tooth. **Before use, adjust the toggle to suit the mouth quadrant in which the Pro-Matrix is being used.**



Fit the band around the tooth and tighten by turning the thumbscrew clockwise. Do not over tighten. For larger teeth, unscrew the band counterclockwise.



Apply wedges and shape the band using a suitable burnishing tool if required.



Perform restoration using normal clinical procedure.



To remove the band, first loosen it by turning the handle 2-3 turns counterclockwise.